



SPECIAL FEATURES

- IP68 rated medical keyboard
- · Nano Silver additive with antimicrobial properties
- Adjustable brightness of backlight
- Fully washable and disinfectable
- 110 keys incl. num pad and track pad

- Waterproof silicone USB protection cap
- Excellent tactile feeling
- ON/OFF switch for keyboard power and backlight
- EN60601-1-2 certified

SPECIFICATION

		Interfaces:	USB
Number of keys:	110 (US-Layout)	Cable length:	2 meters
Keyswitch technology:	membrane	Operating temperature:0℃~ +55℃	
Switching cycles:	approx. 2 mill. (per key)	Storage temperature:	-30°C~+70°C
Protection level:	IP68	Dimensions:	410(L) x 140(W) x 13(T) mm
Housing material:	silicone	Weight:	approx. 1005 g
Housing colour:	white & black	Certification:	CE, FCC, RoHs
Backlight:	blue colour, adjustable brightness	Warranty:	1 year

ORDER INFORMATION

Cat. No.	Product description	Colour	Layout
KSI-U10080	Silicon Touchpad Keyboard	White	US
KSI-U20080	Silicon Touchpad Keyboard	Black	US

Other layouts, configurations and interfaces on request.

GETT Gerätetechnik GmbH

Mittlerer Ring 1 08233 Treuen

Version 04 2017

 Miller of Fung I

 08233 Treuen

 Germany

 Tel:

 +49 37468 660 0

 E-Mail:

 Info@gett-group.com

 Internet:

 www.gett-group.com

GETT North America LLC Suite 203M 2033 S. Voss Road, Houston, TX 77057 United States of America Tel: +832 649 7822 E-Mail: infr@ggett-na.com Internet: www.gett-na.com

 GETT Asia Ltd.

 Unit 8, 5/F, CRE centre,

 No. 889 Cheung Sha Wan Road

 Kowloon, Hong Kong

 Tel:
 +852.3461 9789

 E-Mail:
 info@geta-asia.com

 Internet:
 www.gett-asia.com

gett-asia.com gett-group.com

* EN60601-1-2 CERTIFIED

IEC 60601-1-2 applies to the basic safety and essential performance of Medical Equipment (ME) equipment and ME systems in the presence of electromagnetic disturbances and to electromagnetic disturbances emitted by ME equipment and ME systems.

SHORT OPERATING INSTRUCTIONS



Turn the keyboard on and off by pressing both keys (duration 5 sec).



Turn the backlight on and off with this button.



Adjust the brightness level (down / up) with these buttons.

CLEANING INSTRUCTION

For further details please refer to below link: http://www.gett-group.com/content/service-cleaning-instructions

GETT Gerätetechnik GmbH

2

Version (2017

Mittlerer Ring 1 08233 Treuen Germany Tel: E-Mail: Internet: +49 37468 660 0 info@gett-group.com www.gett-group.com

GETT North America LLC

Suite 203M Suite 203M 2323 S. Voss Road, Houston, TX 77057 United States of America Tel: +832 649 7822 E-Mail: info@gett-na.com Internet: www.gett-na.com

GETT Asia Ltd.

Unit 8, 5/F, CRE centre, No. 889 Cheung Sha Wan Road Kowloon, Hong Kong Tel: +852 3461 9789 E-Mail: info@gett-asia.co info@gett-asia.com www.gett-asia.com Internet:

gett-asia.com gett-group.com